

Invisible treatment  
for visibly beautiful teeth

# Clinical Indications

## for the treatment with Movendi<sup>®</sup>

Spaced teeth	Up to 3 mm between two teeth Up to 4–6 mm for the area of 4–4 in total
Crowding	Up to 2-3 mm with overlapping of two teeth
Rotations	Up to 30° with central incisors and bicuspid Up to 20° with lateral incisors
Intrusion/Extrusion	Up to 1.5 mm per jaw Overbite correction up to 3 mm
Occlusion	No correction possible Overjet up to 3 mm